



Datatrace Time Keeping Certification

Mesa Laboratories, Inc. certifies that all products it designs and sells have been extensively validated and verified to assure that they meet or exceed all design, performance and accuracy specifications. The time keeping accuracy of the circuits used in the Datatrace Tracer were validated during their design and testing to less than 0.05% deviation across the product operating temperature range. Tests were performed using NIST traceable instrumentation. The quartz crystal employed in Datatrace tracers is specified at better than +/- 30ppm initial accuracy. The manufacturing test, and calibration processes and procedures utilized in manufacturing and servicing Datatrace tracers assures that each tracer keeps accurate time.

By  _____ Date 9 June 2006
Ward M. Clark
Manager, Quality & Regulatory Affairs

Mesa Laboratories, Inc.
12100 West 6th Avenue
Lakewood, Colorado 80228
(303) 987-8000
FAX (303) 987-8989